


Product Catalogue

Plastic Surgery



INDEX

SILIMED	04
Certifications.....	05

BREAST IMPLANTS	06
BioDesign Line.....	07
Soft Plus.....	19
Standard Line.....	20

TISSUE EXPANDERS	
BioDesign Expander.....	24
Tissue Expander.....	26

BODY CONTOUR IMPLANTS	28
------------------------------	-----------

FACIAL IMPLANTS	30
------------------------	-----------

MEDGEL	32
---------------	-----------

OUR PROGRAMS	
PSPS.....	34
Silimed Academy.....	35

SILIMED

Silimed's vocation is to contribute to plastic surgery, promoting both research and scientific studies, and developing technologies to meet different needs.



Silimed has been in the market since 1978:

- It is the first company in the world to identify each implant with an individual serial number;
- It is the largest manufacturer of silicone implants in Latin America;
- It is the leader in sales in the Brazilian market;
- It is present in several countries.

Silimed products:

- They have high performance and are developed to better reproduce the anatomy of each body region;
- They are available for the plastic surgery market with the most complete line of silicone implants to meet the demands of various medical specialties.

CERTIFICATES

All Silimed products are produced with quality raw material, which guarantees the necessary certifications to operate in several countries.

ISO 13485

Internationally recognized standard as a quality management system for the medical devices industry.

Brazil

Facilities and production processes inspected and products registered/notified with the National Health Surveillance Agency (ANVISA) and approved, when necessary by the National Institute of Metrology, Quality and Technology (INMETRO). Environmental conditions and production techniques are controlled by ANVISA's Good Manufacturing Practices for Medical Products (GMP).

Europe

Silimed mammary implants have the "CE Mark", which is further proof of the international recognition of the quality of its products.

Other Countries

It acts in compliance with all formal requirements of each Ministry of Health in the countries where its products are distributed.

INFORMATION

- This catalogue was developed to assist plastic surgeons in the choice and knowledge of the Silimed product line.
- The dimensions and volumes described in the tables are presented with approximate values.
- For further information, we recommend that the plastic surgeon read the instructions for use of the product contained in the packaging before using it.
- The silicone may present a color variation ranging from opaque gray to translucent yellow.


MAMMARY IMPLANTS

One of the objectives of Silimed's technological advances is to create implants which offer results which are increasingly suited to patients' needs. Therefore, Silimed offers two lines of mammary implants: BioDesign and Standard, both made up of mechanically resistant silicone elastomer with **low bleed** treatment, which contributes to minimizing rupture and avoiding gel transudation.

SURFACES


PURE POLYURETHANE

- Low capsular contracture rate of 1% in 10 years.¹
- The 3D collagen matrix formed interwoven within the polyurethane foam promotes Velcro-like adhesion.²
- Its adhesive capability reduces the risk of occurring rotation and displacement.²



TRUE TEXTURE

- Own texturing method (*Diffusion Gas Method*) with open pore formation in the shell, without using negative imprinting, sugar or sodium chloride.
- Low capsular contracture rate of 9% in 10 years.³ Other manufacturers report rates of up to 17% in 10 years in primary augmentation.⁴



*Laboratorio de Evaluación y Desarrollo de Biomateriales en el Nordeste de Brasil - Certbio


SMOOTH

- Mechanically resistant silicone elastomer shell with high-performance HSC gel-filled implant.



SOFT PLUS

- It has a smooth and opaque surface, with arithmetic roughness of approximately 5 μm (micrometers), obtained by negative 3D printing.⁵



The BioDesign Collection is an innovative line, developed with exclusive measurements for all shapes and biotypes that offers **5 implant** models with different profiles and **4 different** projections for the same chosen base size. Available* in the following surfaces: Pure Polyurethane and True Texture.

FILLING GEL



- High performance cohesive gel.
- Penetration grade: more consistent
(The lower the grade measured, the more consistent the gel)
- Highest index of gel integration with the shell.⁶
- Maximum gel elasticity index.⁶
- Maintains shape stability.⁶

*The Advance model is only available in Pure Polyurethane surface.

SHAPES



ADVANCE

Round Base • Conical Profile • Medium Pole

- Super projection
- For high chests and young breasts
- For reconstruction after mastectomy surgery



MAXIMUM

Round Base • Spherical Profile • Medium Pole

- Higher projection
- For accentuated breasts
- For reconstruction after mastectomy surgery



NATURAL

Round Base • Tear Drop Profile • Medium Pole

- Natural mammary augmentation and correction
- For reconstruction after mastectomy surgery



ENHANCE

Oval Base • Anatomic Profile • Superior pole

- Upper pole filling
- For tall, thin women with long, narrow thorax
- For reconstruction after mastectomy surgery



NUANCE

Oval Base • Anatomic Profile • Inferior pole

- Lower pole filling
- For women with short and wide thorax
- For reconstruction after mastectomy surgery

Available
in surfaces:


- ◆ Pure Polyurethane
- ◆ True Texture
- ◆ Soft Plus*

*Soft Plus available in round.

PROJECTION

With innovation combined with technology, BioDesign MAXIMUM implants offer greater projection with smaller dimensions and volume.


■ **Greater projection** (Measure C) and **greater arch** length (Measure D) when compared to the competitors' implants with approximate bases. For competitor brands to present the same projection, they need to increase the base dimensions and volume of the implants.



■ Better volume distribution to avoid mechanical complications of an overly heavy implant in the long term.

4 Projections


(LO, MD, HI, XH)
for each base.



PURE POLYURETHANE

REF	VOL	A	B	C	D
	cm ³	cm			
150 XH	150	8,8	8,8	5,0	7,0
140 HI	140			4,5	6,6
180 XH	180	9,3	9,3	5,3	7,4
165 HI	165			4,8	7,0
150 MD	150			4,3	6,6
215 XH	215	9,8	9,8	5,6	7,8
185 HI	185			5,1	7,4
175 MD	175			4,6	7,0
265 XH	265	10,3	10,3	5,9	8,4
235 HI	235			5,4	8,2
210 MD	210			4,9	7,5
190 LO	190			4,4	7,0
300 XH	300	10,8	10,8	6,2	8,8
275 HI	275			5,7	8,5
245 MD	245			5,2	8,0
225 LO	225			4,7	7,4
345 XH	345	11,3	11,3	6,5	9,2
315 HI	315			6,0	8,9
300 MD	300			5,5	8,3
260 LO	260			5,0	7,8
380 XH	380	11,8	11,8	6,8	9,6
355 HI	355			6,3	9,3
325 MD	325			5,8	8,6
295 LO	295			5,3	8,4
435 XH	435	12,3	12,3	7,1	10,2
400 HI	400			6,6	9,6
375 MD	375			6,1	9,2
350 LO	350			5,6	8,8
505 XH	505	12,8	12,8	7,4	10,6
460 HI	460			6,9	10,0
435 MD	435			6,4	9,8
410 LO	410			5,9	9,3

ADVANCE - BIODESIGN



30535

ROUND BASE
CONICAL PROFILE
MEDIAL POLE



EXTRA HIGH
PROJECTION



HIGH
PROJECTION




MODERATE
PROJECTION



LOW
PROJECTION


REF	VOL	A	B	C	D
	cm ³	cm			
155 HI	155	8,1	8,1	4,4	6,9
140 MD	140	8,3	8,3	4,2	6,5
160 MD	160	8,6	8,6	4,4	7,0
180 HI	180	8,7	8,7	4,6	7,2
215 XH	215	9,0	9,0	5,5	8,0
210 HI	210	9,3	9,3	4,8	7,5
250 XH	250			5,7	8,1
225 HI	225	9,6	9,6	4,9	7,7
190 MD	190			4,1	7,1
275 XH	275			5,8	8,3
240 HI	240	9,9	9,9	5,0	7,8
205 MD	205			4,2	7,2
300 XH	300			5,9	8,5
265 HI	265	10,2	10,2	5,1	8,0
225 MD	225			4,3	7,3
320 XH	320			6,0	8,7
280 HI	280			5,2	8,2
240 MD	240	10,5	10,5	4,4	7,4
195 LO	195			3,6	7,1
350 XH	350			6,1	9,0
305 HI	305			5,3	8,4
255 MD	255	10,8	10,8	4,5	7,5
210 LO	210			3,7	7,3
375 XH	375			6,2	9,3
330 HI	330			5,4	8,6
280 MD	280	11,1	11,1	4,6	7,9
235 LO	235			3,8	7,6
410 XH	410			6,3	9,4
355 HI	355	11,4	11,4	5,5	8,9
425 XH	425			6,4	9,6
380 HI	380			5,6	9,0
335 MD	335	11,7	11,7	4,8	8,4
270 LO	270			4,0	7,9
410 HI	410	12,0	12,0	5,7	9,2
485 XH	485			6,6	9,9
435 HI	435			5,8	9,4
390 MD	390	12,3	12,3	5,0	8,8
320 LO	320			4,2	8,3
565 XH	565			6,8	10,3
495 HI	495			6,0	9,8
440 MD	440	12,9	12,9	5,2	9,2
360 LO	360			4,4	8,7
620 XH	620			7,0	10,7
560 HI	560			6,2	10,2
505 MD	505	13,5	13,5	5,4	9,6
420 LO	420			4,6	9,1

MAXIMUM - BIODESIGN



30622

ROUND BASE
SPHERICAL PROFILE
MEDIAL POLE



EXTRA HIGH
PROJECTION



HIGH
PROJECTION




MODERATE
PROJECTION



LOW
PROJECTION

REF	VOL	A	B	C	D
	cm ³	cm			
185 HI	185	9,6	9,6	4,9	6,7
155 MD	155			4,1	5,5
170 MD	170	9,9	9,9	4,2	5,7
245 XH	245	10,2	10,2	5,9	7,5
185 MD	185			4,3	5,9
155 LO	155			3,5	5,3
270 XH	270	10,5	10,5	6,0	7,7
230 HI	230			5,2	6,9
200 MD	200			4,4	6,1
170 LO	170			3,6	5,5
290 XH	290	10,8	10,8	6,1	7,9
250 HI	250			5,3	7,1
220 MD	220			4,5	6,3
180 LO	180			3,7	5,7
315 XH	315			6,2	8,1
270 HI	270	11,1	11,1	5,4	7,3
245 MD	245			4,6	6,5
195 LO	195			3,8	6,0
360 XH	360	11,7	11,7	6,4	8,3
315 HI	315			5,6	7,6
275 MD	275			4,8	6,7
230 LO	230			4,0	6,3
410 XH	410			6,6	8,6
360 HI	360	12,3	12,3	5,8	7,9
315 MD	315			5,0	7,0
270 LO	270			4,2	6,6
470 XH	470	12,9	12,9	6,8	9,0
410 HI	410			6,0	8,2
360 MD	360			5,2	7,4
310 LO	310			4,4	6,9
530 XH	530			7,0	9,3
465 HI	465	13,5	13,5	6,2	8,5
410 MD	410			5,4	7,8
355 LO	355			4,6	7,2
590 XH	590	14,1	14,1	7,2	9,6
525 HI	525			6,4	8,8
465 MD	465			5,6	8,2
400 LO	400			4,8	7,5

NATURAL - BIODESIGN




30637

ROUND BASE
TEAR DROP PROFILE
MEDIAL POLE



EXTRA HIGH
PROJECTION



HIGH
PROJECTION




MODERATE
PROJECTION



LOW
PROJECTION

REF	VOL	A	B	C	D
	cm ³	cm			
220 XH	220	9,4	10,4	5,1	6,6
190 HI	190			4,3	6,1
250 XH	250	9,9	10,9	5,3	6,9
220 HI	220			4,5	6,3
180 MD	180			3,7	5,8
290 XH	290	10,4	11,4	5,5	7,3
250 HI	250			4,7	6,6
210 MD	210			3,9	6,0
335 XH	335	10,9	11,9	5,7	7,6
285 HI	285			4,9	6,9
240 MD	240			4,9	6,2
375 XH	375	11,4	12,4	5,9	7,9
325 HI	325			5,1	7,2
275 MD	275			4,3	6,4
230 LO	230			3,5	6
420 XH	420	11,9	12,9	6,1	8,3
365 HI	365			5,3	7,5
315 MD	315			4,5	6,7
260 LO	260			3,7	6,2
470 XH	470	12,4	13,4	6,3	8,7
410 HI	410			5,5	7,8
355 MD	355			4,7	7
300 LO	300			3,9	6,5
465 HI	465	12,9	13,9	5,7	8,1
405 MD	405			4,9	7,3
345 LO	345			4,1	6,8
515 HI	515	13,4	14,4	5,9	8,4
450 MD	450			5,1	7,7
385 LO	385			4,3	7,1
580 HI	580	13,9	14,9	6,1	8,8
500 MD	500			5,3	8,1
425 LO	425			4,5	7,4
560 MD	560	14,4	15,4	5,5	8,5
480 LO	480			4,7	7,7
620 MD	620	14,9	15,9	5,7	8,9
535 LO	535			4,9	8,1

ENHANCE - BIODESIGN




30674

OVAL BASE
ANATOMICAL PROFILE
SUPERIOR POLE



EXTRA HIGH
PROJECTION



HIGH
PROJECTION




MODERATE
PROJECTION



LOW
PROJECTION


REF	VOL	A	B	C	D
	cm ³	cm			
205 XH	205	10,4	8,6	4,7	6,4
180 HI	180			4,2	5,9
235 XH	235	10,9	9	4,9	6,6
210 HI	210			4,4	6,1
265 XH	265	11,4	9,4	5,1	6,8
235 HI	235			4,6	6,3
305 XH	305	11,9	9,9	5,3	7,0
270 HI	270			4,8	6,6
250 MD	250			4,3	6,3
220 LO	220			3,8	5,8
350 XH	350	12,4	10,3	5,5	7,2
310 HI	310			5,0	6,9
285 MD	285			4,5	6,5
255 LO	255			4,0	6,1
390 XH	390	12,9	10,7	5,7	7,5
345 HI	345			5,2	7,2
320 MD	320			4,7	6,8
285 LO	285			4,2	6,3
435 XH	435	13,4	11,1	5,9	7,8
395 HI	395			5,4	7,5
360 MD	360			4,9	7,0
325 LO	325			4,4	6,6
485 XH	485	13,9	11,5	6,1	8,1
440 HI	440			5,6	7,8
405 MD	405			5,1	7,3
365 LO	365			4,6	6,8
540 XH	540	14,4	11,9	6,3	8,4
490 HI	490			5,8	8,1
450 MD	450			5,3	7,5
410 LO	410			4,8	7,1
595 XH	595	14,9	12,3	6,5	8,7
535 HI	535			6,0	8,4
500 MD	500			5,5	7,8
455 LO	455			5	7,3
555 MD	555	15,4	12,7	5,7	8,2
505 LO	505			5,2	7,6

NUANCE - BIODESIGN




30646

OVAL BASE
ANATOMICAL PROFILE
INFERIOR POLE



EXTRA HIGH
PROJECTION



HIGH
PROJECTION



MODERATE
PROJECTION




LOW
PROJECTION

REF	VOL	A	B	C	D
	cm ³	cm			
255 XH	255			5,6	7,8
225 HI	225	9,2	9,2	4,8	7,2
190 MD	190			4,0	6,7
155 LO	155			3,2	6,2
280 XH	280			5,7	8,0
245 HI	245	9,5	9,5	4,9	7,3
205 MD	205			4,1	6,8
170 LO	170			3,3	6,4
305 XH	305			5,8	8,4
270 HI	270	9,8	9,8	5,0	7,6
225 MD	225			4,2	7,0
185 LO	185			3,4	6,6
325 XH	325			5,9	8,6
285 HI	285	10,1	10,1	5,1	7,8
240 MD	240			4,3	7,2
200 LO	200			3,5	6,8
360 XH	360			6,0	8,8
305 HI	305	10,4	10,4	5,2	8,1
255 MD	255			4,4	7,4
215 LO	215			3,6	7,0
380 XH	380			6,1	9,0
335 HI	335	10,7	10,7	5,3	8,3
280 MD	280			4,5	7,6
240 LO	240			3,7	7,2
410 XH	410			6,2	9,1
365 HI	365	11,0	11,0	5,4	8,5
305 MD	305			4,6	7,8
260 LO	260			3,8	7,4
430 XH	430			6,3	9,2
390 HI	390	11,3	11,3	5,5	8,6
330 MD	330			4,7	8,0
275 LO	275			3,9	7,5
470 XH	470			6,4	9,4
420 HI	420	11,6	11,6	5,6	8,8
350 MD	350			4,8	8,2
305 LO	305			4,0	7,7
500 XH	500			6,5	9,6
445 HI	445	11,9	11,9	5,7	9,0
385 MD	385			4,9	8,4
335 LO	335			4,1	8,0
570 XH	570			6,7	10,0
505 HI	505	12,5	12,5	5,9	9,4
440 MD	440			5,1	8,9
360 LO	360			4,3	8,3
625 XH	625			6,9	10,4
570 HI	570	13,1	13,1	6,1	9,8
505 MD	505			5,3	9,4
430 LO	430			4,5	8,7


TRUE TEXTURE

MAXIMUM - BIODESIGN



20622

ROUND BASE
SPHERICAL PROFILE
MEDIAL POLE



EXTRA HIGH
PROJECTION



HIGH
PROJECTION



MODERATE
PROJECTION




LOW
PROJECTION

CONSULT YOUR SILMED DISTRIBUTOR REGARDING THE DELIVERY
TIME AND AVAILABILITY OF PRODUCTS IN YOUR COUNTRY

MAMMARY IMPLANTS | 15


REF	VOL	A	B	C	D
	cm ³	cm			
250 XH	250	9,8	9,8	5,8	7,3
220 HI	220			5,0	6,6
190 MD	190			4,2	6,0
155 LO	155			3,4	5,3
270 XH	270	10,1	10,1	5,9	7,5
235 HI	235			5,1	6,8
205 MD	205			4,3	6,2
170 LO	170			3,5	5,5
290 XH	290	10,4	10,4	6,0	7,7
255 HI	255			5,2	7,0
220 MD	220			4,4	6,4
185 LO	185			3,6	5,7
320 XH	320	10,7	10,7	6,1	7,9
270 HI	270			5,3	7,2
245 MD	245			4,5	6,6
200 LO	200			3,7	6,0
340 XH	340	10,9	10,9	6,2	8,1
295 HI	295			5,4	7,3
260 MD	260			4,6	6,7
215 LO	215			3,8	6,2
365 XH	365	11,3	11,3	6,3	8,2
320 HI	320			5,5	7,5
275 MD	275			4,7	6,8
230 LO	230			3,9	6,3
380 XH	380	11,6	11,6	6,4	8,4
340 HI	340			5,6	7,7
300 MD	300			4,8	7,0
250 LO	250			4,0	6,5
425 XH	425	11,9	11,9	6,5	8,5
370 HI	370			5,7	7,8
320 MD	320			4,9	7,1
275 LO	275			4,1	6,6
475 XH	475	12,5	12,5	6,7	8,8
425 HI	425			5,9	8,1
360 MD	360			5,1	7,4
315 LO	315			4,3	6,9
550 XH	550	13,1	13,1	6,9	9,1
485 HI	485			6,1	8,4
415 MD	415			5,3	7,7
360 LO	360			4,5	7,2
610 XH	610	13,7	13,7	7,1	9,4
550 HI	550			6,3	8,8
470 MD	470			5,5	8,0
410 LO	410			4,7	7,5

NATURAL - BIODESIGN



20637

ROUND BASE
TEAR DROP PROFILE
MEDIUM POLE



EXTRA HIGH
PROJECTION



HIGH
PROJECTION




MODERATE
PROJECTION



LOW
PROJECTION


REF	VOL	A	B	C	D
	cm ³	cm			
215 XH	215	9,0	10,0	5,0	6,5
190 HI	190			4,2	6,0
250 XH	250	9,5	10,5	5,2	6,8
215 HI	215			4,4	6,2
180 MD	180			3,6	5,4
290 XH	290	10,0	11,0	5,4	7,1
250 HI	250			4,6	6,4
210 MD	210			3,8	5,6
330 XH	330	10,5	11,5	5,6	7,4
285 HI	285			4,8	6,6
240 MD	240			4,0	5,8
375 XH	375	11,0	12,0	5,8	7,7
320 HI	320			5,0	6,9
275 MD	275			4,2	6,1
225 LO	225			3,4	5,6
420 XH	420	11,5	12,5	6,0	8,0
365 HI	365			5,2	7,2
310 MD	310			4,4	6,4
260 LO	260			3,6	5,8
465 XH	465	12,0	13,0	6,2	8,3
410 HI	410			5,4	7,5
350 MD	350			4,6	6,7
295 LO	295			3,8	6,1
520 XH	520	12,5	13,5	6,4	8,6
460 HI	460			5,6	7,8
400 MD	400			4,8	7,0
340 LO	340			4,0	6,4
510 HI	510	13,0	14,0	5,8	8,1
450 MD	450			5,0	7,3
380 LO	380			4,2	6,7
570 HI	570	13,5	14,5	6,0	8,4
500 MD	500			5,2	7,6
425 LO	425			4,4	7,0
635 HI	635	14,0	15,0	6,2	8,7
560 MD	560			5,4	8,0
475 LO	475			4,6	7,3
615 MD	615	14,5	15,5	5,6	8,4
530 LO	530			4,8	7,6

ENHANCE - BIODESIGN




20674

OVAL BASE
ANATOMICAL PROFILE
SUPERIOR POLE



EXTRA HIGH
PROJECTION



HIGH
PROJECTION




MODERATE
PROJECTION



LOW
PROJECTION


REF	VOL	A	B	C	D
	cm ³	cm			
220 XH	220	10,0	8,2	4,6	6,0
185 HI	185			4,1	5,7
240 XH	240	10,5	8,6	4,8	6,3
210 HI	210			4,3	6,0
280 XH	280	11,0	9,0	5,0	6,6
245 HI	245			4,5	6,3
325 XH	325	11,5	9,5	5,2	6,9
285 HI	285			4,7	6,6
250 MD	250			4,2	5,8
220 LO	220			3,7	5,3
350 XH	350	12,0	9,9	5,4	7,1
315 HI	315			4,9	6,9
280 MD	280			4,4	6,1
250 LO	250			3,9	5,9
405 XH	405	12,5	10,3	5,6	7,4
355 HI	355			5,1	7,2
320 MD	320			4,6	6,4
285 LO	285			4,1	5,9
440 XH	440	13,0	10,7	5,8	7,7
395 HI	395			5,3	7,5
360 MD	360			4,8	6,7
320 LO	320			4,3	6,2
485 XH	485	13,5	11,1	6,0	8,0
450 HI	450			5,5	7,8
400 MD	400			5,0	7,0
365 LO	365			4,5	6,5
500 HI	500	14,0	11,5	5,7	8,1
450 MD	450			5,2	7,3
405 LO	405			4,7	6,8
455 LO	455	14,4	11,9	4,9	7,1
545 HI	545	14,5	11,9	5,9	8,4
500 MD	500			5,4	7,6
555 MD	555	15,0	12,3	5,6	7,9
505 LO	505			5,1	7,4

NUANCE - BIODESIGN




20646

OVAL BASE
ANATOMICAL PROFILE
INFERIOR POLE



EXTRA HIGH
PROJECTION



HIGH
PROJECTION



MODERATE
PROJECTION



LOW
PROJECTION

REF	VOL	A	B	C	D
	cm ³	cm			
225 HI	225			4,4	7,2
195 MD	195	9,8	9,8	3,9	6,7
160 LO	160			3,2	6,2
250 HI	250			4,5	7,3
210 MD	210	10	10	4,0	6,8
175 LO	175			3,3	6,4
275 HI	275			4,7	7,6
230 MD	230	10,3	10,3	4,1	7,0
190 LO	190			3,4	6,6
290 HI	290			4,8	7,8
245 MD	245	10,5	10,5	4,2	7,2
210 LO	210			3,5	6,8
315 HI	315			4,9	8,1
265 MD	265	10,9	10,9	4,3	7,4
225 LO	225			3,6	7,0
345 HI	345			5,1	8,3
290 MD	290	11,2	11,2	4,5	7,6
250 LO	250			3,7	7,2
375 HI	375			5,2	8,5
315 MD	315	11,5	11,5	4,6	7,8
265 LO	265			3,9	7,4
400 HI	400			5,3	8,6
340 MD	340	11,8	11,8	4,7	8,0
285 LO	285			4,0	7,5
425 HI	425			5,4	8,8
365 MD	365	12,1	12,1	4,8	8,2
315 LO	315			4,1	7,7
455 HI	455			5,5	9,0
390 MD	390	12,4	12,4	4,9	8,4
340 LO	340			4,2	8,0
520 HI	520			5,6	9,4
450 MD	450	12,9	12,9	5,1	8,9
380 LO	380			4,3	8,3
575 HI	575			5,7	9,8
515 MD	515	13,5	13,5	5,3	9,4
445 LO	445			4,4	8,7
575 MD	575	14,1	14,1	5,4	9,7
495 LO	495			4,5	9,1
580 LO	580	14,7	14,7	5,0	9,8

SOFT⁺Plus


PERFECT COMBINATION OF SOFTNESS AND DESIGN.

SOFT PLUS IS THE RESULT OF THE EVOLUTION OF **THE CLASSIC SMOOTH SURFACE** ADDED TO THE **HSC+ GEL** AND THE **BIODESIGN** CONCEPT.

This new implant has additional benefits compared to smooth implants, such as the complete filling of the breasts.


Soft Plus heralds **the new generation of smooth implants** on the market!

MAXIMUM - BIODESIGN



22622

ROUND BASE
SPHERICAL PROFILE
MEDIAL POLE



HIGH
PROJECTION



MODERATE
PROJECTION



LOW
PROJECTION

STANDARD


This line is a classic option for round mammary implants. Silimed offers two surface options: True Texture and Smooth.

FILLING GEL



- High performance cohesive gel;
- Penetration grade: softer;
(The lower the measured grade, the higher the gel consistency)
- Higher gel integration index with the shell;⁶
- Higher gel fracture toughness index.⁶

SHAPE



MAXIMUM

Round base • Spherical profile • Medium pole

- Mammary correction and augmentation
- Fully filled breasts in a more natural way
- They have an optimum base and projection size for a classic option


Available in
the surface:

- ◆ True Texture
- ◆ Liso

PROJECTION

3 Projections

Extra High (EH), High (H),
and Extra Low (EL)




Line STANDARD

LAUNCH OF A **SILIMED PRODUCT – EXCLUSIVE TO THE EUROPEAN MARKET**

the **TRUE TEXTURE Extra Low Breast Implant (20610)** from Silimed's STANDARD Line, designed to meet the demands of **broader base** and **lower projection**.


REF	VOL	A	B	C
	cm ³	cm		
080	080	8,0	8,0	2,2
100	100	8,6	8,6	2,3
120	120	9,1	9,1	2,4
140	140	10,1	10,1	2,4
160	160	10,5	10,5	2,5
190	190	11,4	11,4	2,6
220	220	12,1	12,1	2,7
250	250	12,3	12,3	2,8
280	280	13,0	13,0	2,9
310	310	13,4	13,4	3,0
340	340	13,7	13,7	3,1
370	370	14,5	14,5	3,1
400	400	14,8	14,8	3,2
450	450	15,4	15,4	3,2
500	500	15,9	15,9	3,3
550	550	16,5	16,5	3,5
600	600	17,0	17,0	3,7
650	650	17,5	17,5	3,8
700	700	17,9	17,9	3,8

MAXIMUM EXTRA LOW



20610


WIDE ROUND BASE
SPHERICAL PROFILE
MEDIUM POLE
EXTRA LOW PROJECTION



Line STANDARD


REF	VOL	A	B	C
	cm ³	cm		
115	115	8,1	8,1	3,2
135	135	8,7	8,7	3,3
155	155	8,8	8,8	3,6
175	175	9,3	9,3	3,7
195	195	9,7	9,7	3,9
215	215	10,0	10,0	3,9
235	235	10,5	10,5	3,9
255	255	11,0	11,0	3,9
285	285	11,4	11,4	4,0
305	305	11,7	11,7	4,0
325	325	11,8	11,8	4,4
355	355	12,1	12,1	4,4
385	385	12,5	12,5	4,5
435	435	13,2	13,2	4,6
485	485	13,6	13,6	4,8
525	525	13,8	13,8	5,0
575	575	14,1	14,1	5,2

MAXIMUM



20621

ROUND BASIS
SPHERICAL PROFILE
MEDIAL POLE
HIGH PROJECTION




REF	VOL	A	B	C
	cm ³	cm		
205 EA	205	9,5	9,5	4,3
225 EA	225	9,7	9,7	4,5
245 EA	245	10,1	10,1	4,5
265 EA	265	10,4	10,4	4,6
285 EA	285	10,8	10,8	4,7
305 EA	305	10,9	10,9	4,8
330 EA	330	11,1	11,1	5,0
355 EA	355	11,5	11,5	5,1
380 EA	380	11,7	11,7	5,4
410 EA	410	12,0	12,0	5,5
450 EA	450	12,0	12,0	5,9
500 EA	500	12,4	12,4	6,3
550 EA	550	12,8	12,8	6,6
600 EA	600	13,1	13,1	6,8
650 EA	650	13,5	13,5	7,0

MAXIMUM

20621 EA


ROUND BASE
SPHERICAL PROFILE
MEDIAL POLE
EXTRA HIGH PROJECTION



Line STANDARD


REF	VOL	A	B	C
	cm ³	cm		
155	155	9,7	9,7	2,9
175	175	9,9	9,9	3,0
195	195	10,1	10,1	3,1
215	215	10,2	10,2	3,4
235	235	10,7	10,7	3,4
255	255	11,0	11,0	3,5
275	275	11,1	11,1	3,7
295	295	11,3	11,3	3,8
315	315	11,4	11,4	4,0
335	335	11,8	11,8	4,1
355	355	11,9	11,9	4,2
385	385	12,3	12,3	4,3
415	415	12,7	12,7	4,4
455	455	12,8	12,8	4,5
495	495	13,3	13,3	4,7

MAXIMUM



10521

ROUND BASE
SPHERICAL PROFILE
MEDIAL POLE
HIGH PROJECTION




REF	VOL	A	B	C
	cm ³	cm		
175 EA	175	8,9	8,9	3,9
195 EA	195	9,2	9,2	4,2
215 EA	215	9,5	9,5	4,4
235 EA	235	9,7	9,7	4,6
255 EA	255	9,8	9,8	4,7
275 EA	275	10,2	10,2	4,8
295 EA	295	10,3	10,3	5,0
315 EA	315	10,5	10,5	5,2
335 EA	335	10,7	10,7	5,3
355 EA	355	10,8	10,8	5,4
385 EA	385	11,0	11,0	5,7
415 EA	415	11,2	11,2	6,0
455 EA	455	11,6	11,6	6,2
495 EA	495	12,1	12,1	6,3
545 EA	545	12,5	12,5	6,5
605 EA	605	12,9	12,9	6,7
655 EA	655	13,1	13,1	6,8
695 EA	695	13,4	13,4	6,9

MAXIMUM

10521 EA

ROUND BASE
SPHERICAL PROFILE
MEDIAL POLE
EXTRA HIGH PROJECTION



TISSUE EXPANDERS

BIODESIGN

The **BIODESIGN tissue expander** line represents the evolution of expanders for breast reconstruction, offering solutions that deliver results increasingly closer to the natural female anatomy.

The shapes of the expanders were specifically designed to complement the various profiles and surfaces of the breast implants in the **BIODESIGN Collection line**.

POLYESTER REINFORCEMENT

The **BIODESIGN tissue expanders** feature an internal polyester reinforcement, designed to maintain the desired shape during the expansion period. This ensures that the expanded area closely matches the corresponding breast implant, while reducing pressure on the tissues and bones located behind the expander base.

REMOTE VALVE

System designed with **2 remote valves**, one adult and one pediatric, externally connected by a tube and a connector pin.

For information on the assembly **sequence of the remote valves**, please refer to our product Instructions for Use (IFU).

REF	BASE		PROJECTION
20800-300-3	10,6	10,6	5,0
20800-400-3	11,2	11,2	5,5
20800-500-3	12,2	12,2	6,2
20800-600-3	13,2	13,2	6,5
20800-700-3	15,0	15,0	6,1



MAXIMUM

TISSUE EXPANDER

20800


ROUND BASE
SPHERICAL PROFILE

Correlating the **BIODESIGN tissue expander** with the BIODESIGN breast implants.

BREAST IMPLANT ADVANCE

Round Base • Conical Profile • Medium Pole

AVAILABLE IN THE SURFACE



BREAST IMPLANT MAXIMUM

Round Base • Spherical Profile • Medium Pole

AVAILABLE IN SURFACES



REF	BASE		PROJECTION
207637-255	10,3	10,3	6,0
207637-380	11,8	11,8	6,7
207637-500	13,0	13,0	7,0
207637-575	13,6	13,6	7,4
207637-640	14,1	14,1	7,5




NATURAL

TISSUE EXPANDER

207637

ROUND BASE
TEAR DROP PROFILE

Correlating the **BIODESIGN tissue expander** with the BIODESIGN breast implants.



BREAST IMPLANT NATURAL

Round Base • Tear Drop Profile • Medium Pole

AVAILABLE IN SURFACES



REF	BASE		PROJECTION
207674-260	10,0	11,0	5,4
207674-345	11,0	12,0	6,0
207674-445	12,0	13,0	6,2
207674-495	12,5	13,5	6,6
207674-550	13,0	14,0	6,8
207674-610	13,5	14,5	7,0
207674-690	14,0	15,0	7,2



ENHANCE

TISSUE EXPANDER

207674

OVAL BASE
ANATOMICAL PROFILE

Correlating the **BIODESIGN tissue expander** with the BIODESIGN breast implants.



BREAST IMPLANT ENHANCE

Oval Base • Anatomical Profile • Superior Pole

AVAILABLE IN SURFACES



TISSUE EXPANDERS

These consist of a mechanically resistant silicone elastomer shell, with a valve connected externally by a tube and a connector pin through which periodic injections of sterile, non-pyrogenic isotonic saline solution will be made until the desired expansion is reached.

INDICATIONS

- Post-mastectomy mammary reconstructions;
- Correction of burn sequels;
- Scar revision in general.

POST-MASTECTOMY DIFFERENTIAL

- Does not affect or alter the color, texture, vascularization and innervation of the tissue that was expanded;^{7,8}
- Allows better planning and follow-up of the reconstruction;
- Ensures a more natural, safe and comfortable result for the woman.⁹


TYPES OF REMOTE VALVES

It has 3 valves: round (adult and child) and luer lock.

For information about the assembly **sequence of the remote valves**, see our product instruction for use (IUP).

REF 800	VOL	A	B
	cm ³	cm	
800 - XXX - 3	200	9,6	4,1
	300	10,6	5,0
	400	11,2	5,5
	500	12,2	6,2
	600	13,2	6,5
	700	15,0	6,1

ROUND




VALVE




REF 810	VOL	A	B	C
	cm ³	cm		
810 - XXX - 3	100	6,0	4,6	5,3
	105	8,7	4,1	3,3
	160	11,0	5,1	3,8
	200	6,0	5,0	7,0
	205	9,5	5,5	4,0
	250	9,6	5,9	5,0
	320	14,1	6,2	4,7
	400	13,6	5,5	5,6
	500	15,8	6,9	5,1
	600	18,0	5,0	7,0

RECTANGULAR



VALVE



BODY CONTOURING IMPLANTS

GLUTEAL IMPLANTS


- Produced with a mechanically resistant elastomer shell;
- Filled with high performance 100% medical grade silicone gel, which has a consistency and hardness designed to resemble the tissues of the target area.¹⁰

GLUTEAL QUARTZO

The Gluteal Quartzo implant is indicated for those who have a hip shape, where the height is higher and the width is smaller.

REF	VOL	A	B	C
	cm ³	cm		
200	200	13,7	9,2	2,8
250	250	14,0	9,5	3,5
300	300	14,3	9,7	3,8
350	350	14,9	10,4	4,2
400	400	15,2	10,9	4,6
450	450	15,4	11,2	4,7
500	500	15,8	11,7	4,8

GLUTEAL QUARTZO



13961

OVAL BASE
HIGH PROJECTION
SMOOTH OPAQUE SURFACE

ROUND GLUTEAL

Ideal for those who have a hip shape, where the height and width are similar.

REF	VOL	A	B
	cm ³	cm	
200	200	10,6	3,7
220	220	10,7	3,8
240	240	10,9	3,9
270	270	11,4	3,9
300	300	11,7	4,0
330	330	12,1	4,1
360	360	12,3	4,3

ROUND GLUTEAL



10910

ROUND BASE

TRANSLUCENT SMOOTH SURFACE

CALF IMPLANT

It is indicated to correct the inequalities or irregularities of the legs.

REF	VOL	A	B	C
	cm ³	cm		
085	085	15,0	4,7	2,4
140	140	19,7	5,4	2,5
180	180	23,1	5,8	2,6

CALF IMPLANT



10410

ASYMMETRICAL BASIS

FACIAL

INDICATIONS

- For mentoplasty, aesthetic or reconstructive surgery;
- Used in association with rhinoplasty or rhytidoplasty surgeries.^{11,12}

CHIN IMPLANT

Rounded shape with perforations for fixation of the connective tissue.

REF	NUM/ REF	A	B	C
		cm		
10 - 035	1	3,5	1,6	0,7
10 - 045	2	4,5	1,7	0,8
10 - 050	3	4,8	1,7	0,8
10 - 055	4	5,3	1,8	0,8
10 - 075	5	7,5	1,9	0,8

CHIN IMPLANT



NUMBER/SIZE

- #1 - SMALL
- #2 - MEDIUM SHORT
- #3 - MEDIUM LONG
- #4 - LARGE
- #5 - CUTTABLE



ANATOMICAL CHIN IMPLANT

- Covers the full lateral expansion of the chin;
- It has a reinforced tab on its internal side for the fixation suture.

REF	NUM/ REF	A	B	C
		cm		
15 - 060	1	6,1	1,4	0,6
15 - 065	2	6,4	1,5	0,7
15 - 080	3	7,9	1,5	0,8

CHIN IMPLANT

ANATOMICAL



NUMBER/SIZE

- #1 - SMALL
- #2 - MEDIUM
- #3 - LARGE



NASAL MODELER

DIFFERENTIAL


- For all age groups and nose types (Caucasian, Negroid and Oriental);
- Provides septum maintenance in vertical position;¹³
- Enables nasal breathing in immediate post-operative period, not obstructing the air flow;¹³
- Easy to clean;¹³
- For single and continuous use for approximately 6 months, and may be replaced by a larger size, as the nose grows.¹³

INDICATIONS

- Prophylaxis of stenosis;
- Scar contractures;
- Septoplasty;
- Primary rhinoplasty surgeries (unilateral cleft lip).

REF	NUM	A	B	C
		cm		
1171 - 019	2	1,90	2,60	0,80
1171 - 020	3	2,00	2,70	0,90
1171 - 021	4	2,10	2,80	1,00
1171 - 022	5	2,20	2,90	1,05
1171 - 023	6	2,30	3,00	1,10
1171 - 024	7	2,40	3,10	1,20
1171 - 025	8	2,50	3,20	1,30
1171 - 026	9	2,60	3,30	1,40
1171 - 027	10	2,70	3,40	1,50
1171 - 028	11	2,80	3,50	1,60
1171 - 029	12	2,90	3,50	1,70

NASAL MODELER



MEDGEL

DIFFERENTIAL

- Made from 100% medical grade silicone gel
- Biocompatible
- Flexible
- Reusable and with soft grip
- No components of animal origin
- No active ingredient



INDICATION

- Helps in the treatment and prevention of keloids, common and hypertrophic scars, recent or old, and associated erythema (redness).¹⁴

SEAL OF AUTHENTICITY

- All Medgel* products receive the Silimed brand engraving during production.

MEDGEL COAT

- For uncovered areas, such as hands, neck and face, that do not have contact with clothing.

REF	A	B	QTD
	cm		
1230 - 015	15	1	1

MEDGEL COAT



*except Medgel Coat

MEDGEL PLATE

■ For wide or scattered scars, e.g., from burns.

REF	A	B	QTD
	cm		
1240 - 010 1	10	10	1

MEDGEL PLATE



MEDGEL STRIP

■ For linear scars, such as abdominoplasty, mammary implants, caesarean sections, heart surgery, etc.

REF	A	B	QTD
	cm		
1240 - 011 1	30	2.0	1
1240 - 011 2	30	2.0	2
1240 - 011H	30	2.0	12

MEDGEL STRIP




MEDGEL DISK

■ For areolar scars, both from mammary reduction and mammary augmentation.

REF	A	B	QTD
	cm		
1240 - 012 2	7	2	2

MEDGEL DISK



FOR MORE INFORMATION

visit our website:
www.silimed.com/en/medgel



OUR PROGRAMS

For Patients

SILIMED PRODUCT REPLACEMENT PROGRAM (PSPS)

Silimed offers a product replacement program in cases of rupture due to manufacturing defect or capsular contracture, Baker grades III or IV.

HOW IT WORKS:

RUPTURE OF IMPLANTS BY MANUFACTURING DEFECT;

Lifetime replacement in cases of proven rupture of mammary implants by manufacturing defect.

CAPSULAR CONTRACTURE;

The replacement will be made in cases of capsular contracture, Baker grades III or IV, for the primary surgery of mammary augmentation or for the first mammary reconstruction surgery, considering the following periods counted from the date of surgery:

- Period of up to 6 years (for surgeries performed between September 1, 2014, and May 3, 2017) and up to 10 years (for surgeries performed as of May 3, 2017) in case of mammary implant with textured surface (TRUE TEXTURE).
- Period of up to 10 years in case of mammary implant with surface coated with polyurethane foam (PURE POLYURETHANE).




Access the QR CODE
above for more
information on the
PSPS

Exclusive for the medical community, offering a platform with scientific content, surgical techniques, and Silimed workshops with plastic surgeons.

BUT WHAT IS THE SILIMED ACADEMY PLATFORM?

A dynamic and versatile virtual environment:

- Exclusive access for plastic surgeons;
- Scientific content, including articles and studies;
- On-demand videos for access anytime, anywhere;
- Webinars with masters of plastic surgery.



Access the QR
CODE above
and **REGISTER!**

Silimed Indústria de Implantes Ltda.
Rua Figueiredo Rocha, 374
Rio de Janeiro - RJ - Brasil - CEP: 21240 - 660

Technical Manager:
Marlos de Oliveira e Souza
Chemist - CRQ (III) RJ: 03230137

New Production Plant
Rodovia Washington Luiz, Sítio 20
área A, Sítio 21 área B, Vila Actura,
Duque de Caxias - Rio de Janeiro - Brazil

Technical Manager:
Alexandre dos Santos Zago
Register CRQ-RJ: 032053861

Registros ANVISA:

Breast Implants: 10102180060 /
10102180066 / 10102180002 / 10102180104
Tissue expanders: 10102180019 / 10102180023
Gluteal implants: 10102180081
Calf implants: 10102180009
Facial implants: 10102180089
Nasal Modeler: 10102180076
Medgel: 10102180037
Medgel Coat: 10102180088

Marca CE: 1434_MDD-083/2020
European Authorized Representative:
SILIMED B.V.

Dirección: Bond Eindhoven, luchthavenweg 81, 5657 EA Eindhoven Room 1.11C floor atrium B
silimed.bv@silimed.com.nl



REFERENCIAS: 1. VÁZQUEZ G, Pellón A. Polyurethane-coated silicone gel breast implants used for 18 years. *Aesthetic Plast Surg.* 2007 Jul-Aug;31(4):330-6. 2. FLEMING, D. Polyurethane foam covered breast implants. *Peters, W. et. al. Biomaterials in plastic surgery: Breast Implants.* Oxford, 2012, P. 96 - 120. 3. Stevens WG, Calobrace MB, Alizadeh K, Zeidler KR, Harrington JL, d'Incelli RC. Ten-year Core Study Data for Sientra's Food and Drug Administration-Approved Round and Shaped Breast Implants with Cohesive Silicone Gel. *Plast Reconstr Surg.* 2018 Apr;141(4S Sientra Shaped and Round Cohesive Gel Implants):75-195. 4. Spear SL, Murphy DK: Allergan Silicone Breast Implant U.S. Core Clinical Study Group. Natrelle round silicone breast implants: Core Study results at 10 years. *Plast Reconstr Surg.* 2014 Jun;133(6):1354-1361. 5. Silimed Indústria de Implantes Ltda. Consolidated Report of Mechanical Tests on Smooth Mammary Implants Ref 22622-XXXLO/MD/HL/XH. Rio de Janeiro, Brasil: Silimed 2021. 6. Kinney BM, Jeffers LLC, Ratliff GE, Carlisle DA. Silicone gel breast implants: science and testing. *Plast Reconstr Surg.* 2014 Jul;134(1 Suppl):475-565. doi: 10.1097/PRS.0000000000000349. 7. Han Y, Zhao J, Tao R, Guo L, Yang H, Zeng W, Song B, Xia W. Repair of Craniomaxillofacial Traumatic Soft Tissue Defects With Tissue Expansion in the Early Stage. *J Craniofac Surg.* 2017 Sep;28(6):1477-1480. doi: 10.1097/SCS.0000000000003852. PMID: 28841593. 8. Motamed S, Niazi F, Atarian S, Motamed A. Post-burn head and neck reconstruction using tissue expanders. *Burns.* 2008 Sep;34(6):878-84. doi: 10.1016/j.burns.2007.11.018. Epub 2008 Apr 2. PMID: 18375071. 9. Casella D, Di Taranto G, Marcasciano M, Lo Torto F, Barellini L, Sordi S, Gaggelli I, Roncella M, Calabrese C, Ribuffo D. Subcutaneous expanders and synthetic mesh for breast reconstruction: Long-term and patient-reported BREAST-Q outcomes of a single-center prospective study. *J Plast Reconstr Aesthet Surg.* 2019 May;72(5):805-812. doi: 10.1016/j.jbpts.2018.12.018. Epub 2018 Dec 16. PMID: 30639155. 10. IUP 006-01-LP-EC: Implante Glúteo - Gel de Silicone - Superfície Lisa. 11. Krishna Patel, MD Kathleyn Brandstetter, MD: Solid Implants in Facial Plastic Surgery: Potential Complications and How to Prevent Them: *Facial Plast Surg* 2016; 32:520-53. 12. Mittelman H, Spencer JR, Chirzanowski DS. Chin region: management of grooves and mandibular hypoplasia with alloplastic implants. *Facial Plast Surg Clin North Am.* 2007 Nov;15(4):445-60. vi. 13. Fonseca JLS. Rinossetoquelloplastia primária nas fissuras labiopalatais unilaterais: resultados preliminares. *Revista Brasileira de Cirurgia Plástica*, 83 (1): 29,1993. 14. IUP 019-01-L-N: Medgel - Gel de Silicone para Tratamento de Cicatrizes / Medgel Coat.